# medinose<sup>+</sup>

## Anti-Allergiegerät

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# English

## Contents

|                                                   |           |
|---------------------------------------------------|-----------|
| <b>1 Introduction</b>                             | <b>12</b> |
| Safety Notes                                      |           |
| <b>2 Interesting Facts about Allergies</b>        | <b>13</b> |
| Allergies – Why are they becoming more frequent?  |           |
| Hay fever – A trifling complaint?                 |           |
| Hay fever – How it arises.                        |           |
| What can help to combat this over-reaction?       |           |
| How phototherapy using the medinose can help you. |           |
| The benefits of your phototherapy with medinose   |           |
| <b>3 Starting Up</b>                              | <b>15</b> |
| Packaging / Scope of Supply                       |           |
| Activation / Replacing the Battery                |           |
| <b>4 Phototherapy with medinose</b>               | <b>16</b> |
| Use of the medinose Appliance                     |           |
| Duration and Frequency of Application             |           |
| Preventive Treatment                              |           |
| Phototherapy – safe and effective                 |           |
| <b>5 General</b>                                  | <b>18</b> |
| Cleaning and Care                                 |           |
| Notes on Disposal                                 |           |
| Technical Data                                    |           |
| <b>6 Warranty</b>                                 | <b>20</b> |

## 1 Introduction

Thank you for your confidence.

Valued customer,

Congratulations on purchasing the medinose. You have bought a quality product from medisana®, one of the world's leading healthcare companies. **Treatment with your medinose appliance is very easy and can be carried out almost anywhere.** The applications lasting only approx. 4 minutes can provide you with a significant improvement within the shortest of times. Please read carefully through the appliance's instructions for use and follow the safety instructions to ensure safe operation and the best possible effect before using it for the first time. Keep these instructions for use in a safe place and include them if the appliance is passed on to a third party.

Thank you for purchasing the medinose. We wish you every success with it.

The medinose appliance from medisana® carries the CE symbol (conformity symbol) confirming the appliance's conformity to "Guideline 93/42/EEC of the Council dated 14 June 1993 governing medicinal products". The corresponding conformity declaration is available at medisana's® premises.

medinose is imported into the EC by:

(authorised EU representative in accordance with the MPG)

Trebs B.V. 6467 JE-Kerkrade, Netherlands

Sales are made through: medisana® AG, Itterpark 7-9, D-40724 Hilden

We reserve the right to make technical modifications for the improvement and further development of the product. © by medisana AG, Germany, Status 11/2004 V1

**CE 0473**

## Safety Notes

- The appliance is licensed for home therapy and must not be used for commercial purposes. It cannot take the place of a visit to the doctor.
- Children may only use the appliance under adult supervision.
- Keep the battery out of the reach of children. The battery is not rechargeable. Do not short-circuit the battery or throw it into a fire.
- Hand it in to your dealer or a local collection point for recycling.
- Only use the appliance for the purpose defined in the instructions for use. The right to claim under the warranty will expire in the event of misuse.
- Only use nasal applicators from medisana®.
- Do not under any circumstances attempt to repair the appliance yourself. Contact your specialist dealer.
- Applications with the medinose are generally perceived as pleasant. Cut short the application immediately if you experience sudden pain or the application becomes unpleasant and consult your doctor.
- The appliance will have limited or no effectiveness in the case of chronic inflammations of the nasal mucosa, deformations of the nasal septum or polyps.

## 2 Interesting Facts about Allergies

### Allergies – Why are they becoming more frequent?

Illnesses caused by allergies such as hay fever and allergic reactions to house dust and animal hair have almost doubled in the last few years. It is not yet medically clear why more and more people are suffering from allergic symptoms. Exaggerated hygiene in childhood may hamper the development of the immune system and favour allergies. It is certain, however, that along with heredity increased environmental pollution plays a significant role in the development of allergies.

### Hay fever – A trifling complaint?

Hay fever, otherwise known as allergic rhinitis, is not actually a fever and does not necessarily have anything to do with hay. It is not caused by a virus but by pollen. Hay fever should always be taken seriously. The allergy can spread to the lower respiratory passages, leading to chronic disorders and finally to asthma. Your doctor can provide you with a diagnosis and the correct individual advice for you.

## GB

### **Hay fever – How it arises.**

It is our immune system's job to fend off and render harmless any contaminants that invade our bodies. In the case of an allergy, the immune system over-reacts to these contaminants. Even the tiniest amounts of allergens, such as flower pollen for example, are enough to induce the typical hay fever symptoms:

- Agonising desire to sneeze combined with violent sneezing attacks,
- Swollen, itchy nasal mucosa,
- Itchy, watery, swollen eyes.

The reason for this is an over-reaction of our defence mechanism. As soon as the flower pollen touches the mucous membranes of the nose and eyes, the body mobilises all its defences. Too many antibodies are created. These sit on the so-called mast cells. The mast cells in turn release tissue substances (primarily histamine) that lead to the typical miseries of hay fever.

### **What can help to combat this over-reaction?**

Of course, the easiest way of preventing an allergy is to consistently avoid all substances that trigger allergies. However, in practice this is hardly possible as one can come across allergens not only outside but also at home (due to house dust) or in foodstuffs.

In addition to drug treatment using antihistamines and cortisone, inoculations and acupuncture, specialist centres successfully use cost-intensive phototherapy. In this case red light rays of one wavelength (approx. 660 nm) are generated by special light-emitting diodes to eliminate the allergic symptoms and the over-sensitivity.

### **How phototherapy using the medinose can help you.**

Phototherapy is a new type of allergy treatment without drugs.

Now with your medinose from medisana® you can quite easily carry out yourself what until recently was only available as a very cost-intensive therapy in special centres. Diodes that are able to emit red light of exactly the right wavelength (approx. 660 nm) have been developed for the medinose. Alleviation of allergic symptoms and lessening of the over-sensitivity of the nasal mucosa sets in at this wavelength.

## The benefits of your phototherapy with medinose:

- Phototherapy is a drug-free treatment.
- The highly effective application lasting only approx. 4 minutes can be carried out practically anywhere.
- The treatment has no effect on other organs or body parts due to its use locally.

The typical symptoms of hay fever were significantly reduced in clinical tests using the medinose. Some symptoms were alleviated or even eliminated completely in a matter of only a few hours. In the majority of patients it was possible to identify a clear improvement in the condition within two weeks.

## 3 Starting Up

### Packaging / Scope of Supply

We would like our high-quality medical products to reach you in the same condition as they left us. For this reason your medinose is packed in a sturdy protective package made from environmentally compatible materials. If you do not want to use the box for storing the appliance, it may be sent for recycling with your waste paper. Please dispose of the plastic films at the local waste site. The package contains the following parts:

- 1 medinose appliance
- 1 9 volt monobloc battery 6LR61
- 1 Nasal applicator
- 1 Instructions for Use

Please consult your dealer if you discover that the appliance is damaged or fails to function as expected after unpacking.

### Activation / Replacing the Battery

**Important:** The 9V monobloc battery is already inserted on delivery of the appliance. The contacts are protected with an insulating strip to prevent the battery from discharging. This insulating strip must be removed prior to starting the appliance up for the first time. Press down on the symbol on the battery compartment lid (back of the appliance) and at the same time slide it in the direction of the arrow. Now pull out the insulating strip by the tab. Then replace the battery compartment lid and slide it upwards until it clicks audibly into place.

GB

After approx. 100 applications, the green pilot LED will flash after pressing the ON button if the battery power is insufficient for another application. Open the battery compartment (as described above), remove the spent battery and insert a new 9 volt monobloc battery 6LR61 (in the position shown inside the battery compartment). Close the battery compartment again.

**Caution!** Batteries may not be disposed of with your domestic waste under any circumstances. Hand them in to your dealer or to your local collection point for recycling.

## 4 Phototherapy with medinose

### Use of the medinose Appliance

- Insert the plug on the cable of the nasal applicator into the socket of your medinose appliance.
- The nasal applicator must be spread slightly apart to allow the LEDs to be inserted in the nostrils. Insert the LEDs as far as possible into the nostrils. The nasal applicator should nevertheless sit comfortably.
- Press the ON button to start the treatment and begin light emission. The green pilot LED on the front of your medinose will be lit constantly during use.
- Light emission will be ended automatically after approx. 4 1/2 minutes. The green pilot LED will go out thus indicating that the application has ended.
- Now remove the LEDs from the nostrils.

### Notes:

- If the green pilot LED flashes this indicates that the cable is not plugged in correctly or that there is insufficient voltage in the battery for an application.
- If the green pilot LED neither lights up nor flashes this indicates that the battery has been incorrectly inserted or is completely discharged.

### Duration and Frequency of Application

- An application with the medinose appliance only lasts approx. 4 minutes. During this time, the LEDs radiate red light having a wavelength of approx. 660 nm and a defined intensity. The automatic shut-off ensures that the duration of treatment always remains the same and prevents the nasal mucosa from being overloaded.

- The medinose should be used three times a day at the beginning of treatment. Two applications a day will be sufficient when the symptoms ease. Treatment may be reduced to once a day if there is continued improvement.
- Treatment should be spread over the course of a day. The first application should be carried out either before or after breakfast and the last in the evening before going to bed whilst the third application can be carried out during the day simply by taking the appliance with you and using it when the opportunity arises.

**Note:**

The applications should be continued for at least one week and preferably two to achieve the desired treatment result.

**Preventive Treatment**

Your medinose appliance may be used not only in the acute phase of the illness but also as a preventive measure. You will most likely know the period during which your allergy starts or the type of pollen that causes it.

Find out about the current pollen count situation, e.g. on the Internet. This way you can gain experience of when the different types of pollen are likely to be most active in your area and take preventive measures in good time.

Begin prophylactic treatment a good week before the pollen that triggers your allergy is forecast. Even at this stage you should carry out an application two to three times a day depending on the severity of your allergy.

**Phototherapy – safe and effective**

The red light used for phototherapy with medinose has no negative effects on the nasal mucosa. Even with therapy continued over a longish period, there is no fear of any structural impairment of the tissue when no more than three applications per day are administered. If your doctor suggests medical treatment due to the severity of your allergy, it may be combined with phototherapy without hesitation. There are no known interactions with anti-allergenic drugs such as antihistamines or corticosteroids.

## 5 General

### Cleaning and Care

- Clean both LEDs after each use with a damp cloth and a mild soap.
- Do not use aggressive detergents or stiff brushes under any circumstances.
- Do not immerse either the appliance or the cable with the LEDs in water. If liquid has managed to get into the appliance, remove the battery immediately and contact your specialist dealer.
- The nasal applicator may be disinfected with an alcoholic disinfectant, e.g. 70% ethanol.
- You may purchase additional nasal applicators from your specialist dealer so that your medinose appliance can be used by several people under hygienic conditions.
- Remove the battery from the appliance if you do not intend to use it for a longish period. This will prevent the battery from leaking.
- No further maintenance or checking is necessary if the appliance is used in accordance with these instructions for use.

### Notes on Disposal

- The packaging consists entirely of environmentally compatible materials that should be disposed of at the local waste site as valuable secondary raw materials. The cardboard box may be disposed of with your waste paper. Please dispose of the plastic films at the local waste site.
- Dispose of the appliance in an environmentally friendly manner in accordance with the statutory regulations when it reaches the end of its service life. Remove the battery first and hand it in to your dealer or the local collection point for recycling.

## Technical Data

|                           |                                                                                       |
|---------------------------|---------------------------------------------------------------------------------------|
| Model:                    | medinose                                                                              |
| Application:              | Phototherapy for allergic reactions of the nasal mucosa                               |
| Duration of application:  | 4 1/2 minutes +/- 10%<br>Automatic shut-off after recommended duration of application |
| Frequency of application: | Maximum 3 applications per day                                                        |
| Light emission:           | Red light of approx. 660 nm at<br>low intensity of 6 mW per diode                     |
| Storage conditions:       | -20 °C to +50 °C, maximum 85% relative humidity,<br>non-condensing                    |
| Operating conditions:     | +10 °C to +40 °C, maximum 70% relative humidity,<br>non-condensing                    |
| Dimensions:               | approx. 90 x 65 x 30 mm                                                               |
| Weight:                   | approx. 120 g including battery                                                       |
| Battery:                  | 9 volt type 6LR61                                                                     |
| Service life of battery:  | approx. 100 applications                                                              |
| Appliance classification: | BF                                                                                    |
| Article No.               | 45020                                                                                 |

**CE 0473**

## 6 Warranty

### Warranty and Repair Conditions

In the case of some faults or uncertainties it is helpful to speak to a specialist. We would ask you, therefore, to consult your specialist dealer first of all in the event of a warranty claim. However, if it is necessary despite this to send the appliance back to medisana®, please make a note of the fault and enclose it with the warranty card and/or the receipt.

### The following warranty conditions are applicable in this case:

1. medisana® provides a 3-year replacement warranty on products from the date of purchase. Evidence of the date of purchase in the event of a warranty claim must be provided by way of the properly filled out warranty card or the receipt.
2. Defects resulting from material or manufacturing faults will be remedied free of charge within the warranty period.
3. A service rendered under the terms of the warranty will not result in any extension of the warranty period for either the appliance or any parts replaced.
4. The following are excluded from the warranty:
  - a) All damage arising as a result of improper treatment, e.g. not observing the instructions for use.
  - b) Damage that is attributable to repairs or intervention by the purchaser or unauthorised third parties.
  - c) Transport damage that has arisen en route from the manufacturer to the consumer or when the appliance has been sent to the customer service centre.
  - d) Accessories that are subject to normal wear and tear (battery, etc.).
5. Liability for direct or indirect consequential damage that is caused by the appliance is also excluded if the damage to the appliance is acknowledged as being a warranty claim.

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